

## Pocket-Tech PT-150 and PT-300



Order # **116-PT-150**

Order # **116-PT-300**

### Product Information:

Model: **PT-150**  
Series: Pocket-Tech Balances  
Capacity: 150 grams  
Readability: 0.1 g  
Linearity: +/- 0.1g  
Weighing Modes: g, oz, dwt, ozt  
Tare Range: To capacity by subtraction  
Calibration: Auto calibration through keypad and 100 gram weight  
Stabilization Time (sec): 3 seconds  
Maximum Overload: 150% of capacity  
Overload Indication:  
Underload Indication: Negative Value  
Auto Off: 4 minutes of inactivity  
Memory Function: Allows cumulative sample weighing  
Optimum Temperature: 64-77 F (18-25 C)  
Optimum Humidity: Less than 80% relative humidity  
Display Type: 0.50  
Low Battery Indication: yes  
Power Requirements: 3 x AAA Alkaline Batteries  
Weighing Platform Size: 2.75 x 3.1 inches  
Dimensions (hwxwd): 0.8 x 5.25 x 1.25 inches  
Net Weight: 177 grams / 6.2 oz  
Warranty: **New 5YEAR WARRANTY**  
Shipping Weight: 1 lb./455 gram  
Box Dimensions: 8 x 5.25 x 1.25 inches

The **Pocket-Tech** series have a built-in heavy-duty hinged hard case with four weighing modes: grams, ounces, troy ounces and pennyweights.

It has digital auto calibration, advanced memory function and turns on to last weighing mode.

It has 4 dye easy operation and battery saving auto-off. This balance does not have an adapter port, so is strictly for battery operation.

Five Year Warranty. A versatile pocket balance.

### Product Information:

Model: **PT-300**  
Series: Pocket-Tech  
Capacity: 300g  
Readability: 0.1g  
Linearity: +/-0.1g  
Weighing Modes: g,oz,dwt,ozt  
Tare Range: to capacity by subtraction  
Calibration: auto calibration through keypad + 300gm wt. Incl.  
Stabilization Time (sec): within 3 seconds  
Maximum Overload: 150% of capacity  
Overload Indication:  
Underload Indication: Negative value  
Auto Off: 4 Minutes of no activity  
Memory Function: Allows cumulative sample weighing  
Optimum Temperature: 64-77F (18-25c)  
Optimum Humidity: < 80% Relative humidity  
Display Type: 0.50  
Low Battery Indication: Yes  
Power Requirements: 3 x AAA Alkaline Batteries  
Weighing Platform Size: 2.75 x 3.1 inches  
Dimensions (hwxwd): 4.9 x 3.3 x 0.8  
Net Weight: 177g/6.2oz  
Warranty: **New FIVE YEAR WARRANTY**  
Shipping Weight: 1.0 lbs./ 455g  
Box Dimensions: 8 x 5.25 x 1.25 inches