

## Acculab Pocket Pro Series 116-PP-2060D



**Description:** Fits in your pocket for use in the lab or the field. This versatile balance has 4 weighing modes, auto keypad calibration with weights included, auto off, and overload and low battery indications.

### Product Information:

Product Type: Acculab  
Model: PP-2060D  
Series: Pocket Pro Series  
Capacity: 60 g x 0.01 g, fine range 20 g x 0.001 g.  
Readability: 0.01 or 0.001 w/in. fine; range  
Linearity: +/- 0.01; or +/- 0.002 w/in Fine; range  
Weighing Modes: g, oz., ct. dwt.  
Tare Range: To capacity by subtraction  
Calibration: Auto from keypad w/included 50-gram weight  
Stabilization Time (sec): 3  
Maximum Overload: 150% of capacity  
Overload Indication: H reading on LCD  
Underload Indication: L reading on LCD  
Auto Off: 4 minutes of inactivity (Battery Only)  
Optimum Temperature: 64-77F (18-25C)  
Optimum Humidity: less than 80% Rh  
Display Type: LCD 40 mm x 10 mm (w x h)  
Low Battery Indication: LO on LCD  
Power Requirements: AC adapter or 9 volt alkaline battery  
Weighing Platform Size: 82 x 62mm  
Dimensions (hwxwd): 35 x90 x 150mm  
Warranty: 2 year  
Other Features: Portable, 2-year warranty