

Acculab SV-Series Bench Scale 116-SV-30,50,100



Each SV-Series bench scale features two weighing modes, kg and lb., parts counting, (5, 10, 15, 20 pieces), autocalibration from keypad with a large stainless steel platform. The remote display head and large angled LCD allow for easy viewing and operation.

Product Information:

Product Type: **Acculab**
 Model: **Please see Specifications below.**
 Series: **Bench Scale**
 Capacity: **Please see Specifications below.**
 Readability: **Please see Specifications below.**
 Linearity: **Please see Specifications below.**
 Weighing Modes: **kg:lb: parts counting presets of 5,10,25 & 50 pcs**
 Tare Range: **To capacity by subtraction**
 Calibration: **Auto calibration from keypad with optional calibration weight**
 Stabilization Time (sec): **3 seconds**
 Maximum Overload: **200.0%**
 Overload Indication: **E with audio tone**
 Underload Indication: **Flashing negative values**
 Auto Off: **4 minutes no activity (battery only)**
 Optimum Temperature: **64 degrees -77 degrees F(18-25c)**
 Optimum Humidity: **less than 80% RH**
 Display Type: **LCD/.75**
 Low Battery Indication: **LO with 10 beep signal**
 Power Requirements: **9 v Alkaline or nicad battery, or AC adapter (incl.)**
 Weighing Platform Size: **12 x 12.5 inch stainless steel platform (305 x 315mm)**
 Dimensions (hxwx): **4 x 14.5 x 18.5 inches**
 Warranty: **2 year**
 Other Features:
 Optional Calibration mass 30 kg, NBS Class F part #116-8006F. Shipping weight for this balance is 15 lbs. and the box dimensions are 22"L x 14.5" w x 4" height.

Specifications

	#116-SV-30	#116-SV-50	#116-SV-100
Capacity	30kg/66lbs.	50kg/110lbs.	100kg/220lbs.
Readability	0.005kg	0.01kg	0.02kg
Linearity	+/- 0.2mg	+/- 0.2mg	+/-0.2mg
Cal. Weight	30kg(Option) #116-8006F	50kg(Option) #116-8006F	100kg(Option) #116-8006F