

Model 400-201 Type A

New!
Shorter Barrel!
Improved Grip!
Precision Ground Spring!

Pencil Style Durometer

Type A: For soft rubber and plastics

The durometer's 1/2" diameter base is the smallest base permitted by ASTM D2240 specifications. This allows it to be used in confined or hard to reach areas other durometers cannot test. These durometers meet or exceeds all aspects of ASTM D2240. The instrument uses a dial gauge with a full 360-degree scale and bold, easy to read numbers. The dial scale can easily be read to 0.5 point accuracy. A standard max-hold hand holds the peak reading until reset. The dial has a red zone to indicate unreliable readings. Readings in this area should not be reported for either Type A or D.

PTC Metrology® is an **A2LA Accredited Laboratory** for Durometer calibration. Certificate number 1896.01



**Pencil-style
0.5"Ø Base**

- Accurate, fast, and affordable testing of soft rubber and plastic.
- Pencil Style
- See Also [306L](#)
- See Also [408](#)
- See Also [211](#)