

e2000™
MILLENNIUM SERIES DUROMETER

Model 400-502D Type D
(US patent 436874)

Type D: For hard rubber and plastics

Specifications

1. Meets or exceeds current ASTM-D 2240 Specifications.
2. Range (ASTM Type D) 0 to 100.
3. Accuracy ± 1 point.
4. Comes Complete with Test Block and Carrying Case
5. Rubber Grip for improved handling stability.



e2000™ MILLENNIUM SERIES DUROMETERS are available in both analog and digital models.

- For hard rubber, plastics, epoxies, soft metals, etc.
- Most Popular D style
- Accurate, fast, and affordable testing of hard rubber and plastic
- See Also 409

Analog Durometers 500 Series

Model	ASTM Type
501A	Type A
501B	Type B
502C	Type C
502D	Type D
501/O	Type O
502DO	Type DO

Digital Durometers 500 Series

Model	ASTM Type
511A	Type A
511B	Type B
512C	Type C
512D	Type D
511/O	Type O
512DO	Type DO