

137-SV-300/500

High Capacity Precision Hand Wheel Test Stand



SV-300 (capacity 300 lbs.)

SV-500 (capacity 500 lbs.)

- **Hand wheel design provides a high degree of precise control**
- **One stand for both compression (push) and tension (pull) tests**
- **Gauge remains in place when hand wheel is released**
- **Stroke 15"**

Force gauge sold separately

Models SV-300/500 have greater clearance to accommodate oversize materials and a long stroke to test materials that require a greater range of movement. SV-300/500 are precise push/pull test stands that are ideal for force applications requiring relatively consistent speed and maximum gauge control. The force gauge remains in position when the wheel is released which allows the user to easily change samples between tests. The hand wheel operation will provide repeatable test results from sample to sample. The SV-300/500 will accommodate any of the standard range Imada force gauges.