

## Test Stands

### Model 400- 473

Deadweight Test Stand with both 1kg & 5kg weights  
Designed for  
ASTM Type A, B, C, D, DO  
and O Durometers

### Model 400- 470

Deadweight Test Stand with 400g weight for Type OO,  
OOO and SL Durometers.  
Designed for Models 302SL, 411 & 412

### Model 400- 471

Deadweight Test Stand with 1kg weight for Type A, B  
and O Durometers.

Designed for Models 306L, 306BL, 408, 408B, &  
410

### Model 400- 475

Deadweight Test Stand with 5kg weight  
For ASTM Type C, D and DO Durometers.

Designed for Models 307L, 307CL, 409, 409C &  
413 Durometers



Design Model for 473,  
470,471 and 475.

## **Model 400-320 Spring Loaded Test Stand**

Designed for A.S.T.M. Type A, B, C, D, DO and O Durometers.



## **Model 400- 441**

Constant Load Weight for Type A, B, or O Durometers

*Accurate Repetitive hand-held testing with durometer Models 306L, 306BL, 408, 408B and 410.*

