

## THE MASTERPIECE COLLECTION OF MAGNIFIERS FEATURE THESE ADVANTAGES

- Optical Glass Lenses ground and polished to precision standards.
- Scientifically balanced and contoured handles for comfortable use in either hand.
- Functional Design with wide rim to protect lens surfaces.
- Individually packaged - makes ideal gifts



### **106-M-800 DP Masterpiece Display**

Display contains: 2 each of M-802, M-803, M-804 and 3 each of M-824 & M-810 magnifiers.(shown above)



### **106-M-810 Illuminating Hand Magnifier**

2"x4", 4.5 diopter lens is ground and polished to precision ophthalmic standards. Bright 6 watt bulb and brilliant reflector cover adapts for both right and left hand use. Natural hold switch activates light when in use and allows for automatic shut-off of light when not in use. Detachable cord with strain relief easily unplugs for use with or without light. Additional bulbs and power supplies can be ordered separately for convenient replacement.(shown above)

### **106-M-804 Round Magnifier**

4" round lens, 4 diopter power, 10" focal length.



### **106-M-824 Rectangular Magnifier**

2"x4" rectangular lens, 4.5 diopter power, 9" focal length.

